TECHNICAL MANUAL

OPERATOR'S, UNIT AND
DIRECT SUPPORT MAINTENANCE MANUAL
(INCLUDING REPAIR PARTS AND SPECIAL TOOLS LIST)
FOR

KOEHLER ELECTRIC PENSKY-MARTENS TESTER

NSN 6630-00-530-0987

PRECISION PENSKY-MARTENS FLASH TESTER

MODEL 74537

NSN 6630-00-244-9415

This technical manual is an authentication of the manufacturer's commercial literature and does not conform with the format and the content requirements normally associated with Army technical manuals. This technical manual does, however, contain all essential information required to operate and maintain the equipment.

Approved for public release; distribution is unlimited.

HEADQUARTERS, DEPARTMENT OF THE ARMY

28 SEPTEMBER 1990

This technical manual is an authentication of the manufacturer's commercial literature and does not conform with the format and the content requirements normally associated with Army technical manuals. This technical manual does, however, contain all essential information required to operate and maintain the equipment.

Approved for public release; distribution is unlimited.

SUPPLEMENTARY INTRODUCTORY MATERIAL

1-1. Maintenance Forms and Records.

Department of the Army forms and procedures used for equipment maintenance will be those described by DA Pam 738-750, The Army Maintenance Management System.

1-2. Reporting Errors and Recommending Improvements.

You can help improve this manual. If you find any mistakes or if you know of a way to improve the procedures, please let us know. Mail your letters, DA Form 2028 (Recommended Changes to Publications and Blank Forms), or DA Form 2028-2 located in the back of this manual, directly to: Commander, U.S. Army Troop Support Command, ATTN: AMSTR-MCTS, 4300 Goodfellow Blvd., St. Louis, MO 63120-1798. A reply will be furnished to you.

1-3. Destruction of Army Material to Prevent Enemy Use.

Refer to TM 750-244-3 for instructions covering the destruction of Army Material to prevent enemy use.

- 1-4. Administrative Storage of Equipment.
- a. Placement of equipment in administrative storage should be for short periods of time when a shortage of maintenance effort exists. Items should be in mission readiness within 24 hours or within the time factors as determined by the directing authority. During the storage period appropriate maintenance records will be kept.
- b. Before placing equipment in administrative storage, current preventive maintenance checks and services should be completed. Shortcomings and deficiencies should be corrected, and all modification work orders (MWO's) should be applied.
- c. Storage site selection. Inside storage is preferred for items selected for administrative storage. If inside storage is not available, trucks, vans, conex containers and other containers may be used.



K16200 & K16270

ELECTRIC PENSKY-MARTENS TESTER

ASTM D93 & E134

FLASH POINT

BY

PENSKY-MARTENS CLOSED TESTER

KOEHLER

K16200 & K16270

ELECTRIC PENSKY MARTENS TESTER

TABLE OF CONTENTS

SECTION:

- (A) ASSEMBLY PROCEDURE
- (B) INSTRUCTIONS FOR USE & SPARE PARTS
- (C) WIRING DIAGRAM
- (D) MECHANICAL ASSEMBLY DRAWING

<u>K16200 & K16270</u> <u>ELECTRIC PENSKY MARTENS TESTER</u>

SAFETY AND HAZARD WARNING

THIS EQUIPMENT MAY INVOLVE HAZARDOUS MATERIAL

AND OPERATIONS. THIS MANUAL DOES NOT PURPORT
TO ADDRESS ALL OF THE SAFETY PROBLEMS ASSOCIATED

WITH THE USE OF THE EQUIPMENT. IT IS THE RESPONSIBILITY

OF WHOEVER USES THIS EQUIPMENT TO CONSULT AND

ESTABLISH APPROPRIATE SAFETY AND HEALTH PRACTICES,

AND DETERMINE THE APPLICABILITY OF REGULATORY

LIMITATIONS PRIOR TO USE.

NOTE: AS A SAFETY PRECAUTION, NEVER
USE UNREGULATED GAS WITH THIS TESTER

SECTION A

(1) ASSEMBLY PROCEDURE:

NOTE: Unit is partially disassembled to fit into preformed carton for shipping.

TO REASSEMBLE: (See diagram K16200 Electric P.M. Tester)
Components numbered on Diagram K16010, K16020. (-3 -11),

(-3 -6) & (-8 -1) have been disassembled for shipping.

- (A) Place the mechanical top (K16010) into the cup (K16020) and place on Air Bath.
- (B) Install the main gas pipe (-3 -6) into the 1/4 pipe tee using a 3/4 Hex open end or adjustable wrench only, (DO NOT USE PIPE WRENCH) and tighten securely.
- (C) Place the cup holder (-3 -11) on the main gas pipe (-3 -6) at a convenient height and tighten thumb screw.
- (D) Install the support rod (-8 1) into the base and tighten set screw in base (-0 1) to secure rod.
 - (E) Install the stirrer motor on the support rod (-8 -1).
- (F) Connect the flex cable from the mechanical top to the motor and adjust the motor to the angle shown in diagram.
- (G) Connect gas inlet to any regulated low pressure (5-10 p.s.i.) gas supply (such as L.P.G. or natural gas). Do not use direct unregulated pressure from an L.P.G. tank.

 Unit is now ready for operation. Proceed to instructions page (1).

ELECTRIC PENSKY MARTENS TESTER (K16200)

(ASTM D93 and API 510)

INSTRUCTIONS FOR USE

- 1. Thoroughly clean and dry all parts of the cup and its accessories before starting the test. Particular care should be taken to avoid the presence of any gasoline or naptha used to clean the apparatus after a previous test.
- 2. Fill the cup with the oil to be tested up to the level indicated by the filling mark.
- 3. Place the lid on the cup and set the latter in the stove. Take care to have the locating devices properly engaged.
- 4. Insert the thermometer. If it is known that the oil will flash above 220°F., the "P.M. High" thermometer may be selected; otherwise it is preferable to start with the "P.M. Low" thermometer, and then change in case a temperature of 220 to 230°F. is reached.
- 5. Light the test flame and adjust by means of the valve screw on the burner block, so that it is 5/32" diameter, the same size as the bead provided for comparison.
- 6. Connect the instrument to the proper source of electric current, with the electrical cord provided.
- 7. Adjust the supply of heat by adjusting the dial on the powerstat until the temperature reading on the thermometer increases by not less than 9°, nor more than 11°F., per minute.
- 8. Connect the stirrer to the stirrer motor. <u>NOTE</u>: Use receptacle in base of powertrol heater and place switch to ON.
- 9. Apply the test flame at each temperature reading, which is a multiple of 2°F., Up to 220°F. For the temperature range above 220°F., apply the test flame at each temperature reading which is a multiple of 5°F. Apply the test flame by operating the knurled hand knob controlling the shutter and test flame burner, so that the flame is lowered in one-half second, left in its lowered position for one second, and quickly raised to its high position. Discontinue stirring during the application of the test flame.
- 10. Record as the flash point the temperature read on the thermometer at the time of the flame application that causes a distinct flash in the interior of the cup. The true flash must not be confused with the bluish halo that sometimes surrounds the test flame for the applications preceding the one that causes the actual flash.

ELECTRIC PENSKY MARTENS TESTER

INSTRUCTIONS FOR USE

(2 of 2)

Samples of asphalts or very viscous materials may be warmed until they are reasonbly fluid before they are tested. However, no sample should be heated more than is absolutely necessary. It shall never be heated above a temperature of 30°F. below its excepted flash point.


For the determination of flash point of suspensions of solids, bring the material to be tested and the tester to a temperature of 60° +/- 10° F., or 20° F. lower than the estimated flash point, whichever is lower. Turn the stirrer 250 +/- 10 RPM, stirring in a downward direction. Raise the temperature throughout the duration of the test at a rate of not less than 2° nor more than 3° F. With the exception of these requirements for rates of stirring and heating, proceed as prescribed above.


For further details consult ASTM D93.


PARTS LIST

K16010	TOP
K16020	CUP
K160-3-3	Brass Bell (K160-3-3)
K160-3-2	Cast Iron Bell (K160-3-2)
K4160-1-14	Thermometer Ferrule Adapter
K145-8	Thermometer Ferrule (Includes K145-8-3 Alum. Ring)
K160-9	Flexible Shaft & Coupling
K160-3-4	Upright Rod (K160-3-4)
K162-2MO	Motor, Stirrer
	Powerstat (280-115-001-115V) (280-230-001-230V)
	Heater, 750W. (225-115-001-115V) or (225-230-001-230V)
K420-0-2	Heater Shell (K420-0-2)
K162-0-1	Base (K162-0-1)

KOEHLER INSTRUMENT COMPANY, INC. Bohemia, L.I. New York









WARRANTY POLICY

Any product* manufactured by Koehler Instrument Co. , Inc. (hereinafter referred to as the company) is sold on the following basis and none other. <u>ALL IMPLIED WARRANTIES OF MERCHANTABILITY AND OF FITNESS FOR A PARTICULAR PURPOSE ARE HEREBY EXPRESSLY EXCLUDED</u>.

The following warranty shall apply, and no other warranty, express or implied, shall apply.

If within one year from date of shipment the product fails because of defective material or poor workmanship, the company will repair or replace, without charge, any product that has failed provided:

- a) the product has been properly installed, operated and maintained.
- b) the company is advised in writing of the malfunction and authorizes the return of the product to the factory.
- c) All transportation charges for the return to the factory are prepaid. (Products will be returned freight collect.)
- d) A complete description of the reason for return must accompany the unit.

NOTE: A nominal handling charge for inspection will be made on units for which a claimed defect cannot be confirmed.


THE COMPANY'S SOLE LIABILITY HEREUNDER SHALL BE TO REPAIR OR REPLACE ANY PRODUCT WHICH HAS NOT COMPLIED WITH THIS WARRANTY.

In no event shall the company be liable for:

- Prospective profits or special, indirect or consequential damages caused by failure of its product.
- 2) Any charges for labor or materials for work done on its products by others.

*Wherever used in this Warranty Policy the term "product" shall mean any items manufactured and/or sold by Koehler Instrument Co., Inc.

Precision[®] Instruction Manual TS-74537 AP-9



Precision Scientific

Precision™ Instruction Manual

TS-74537 AP -9

Pensky-Martens Flash Testers Catalog 74537, 74540 and 74545

Introduction

Your satisfaction and safety are important to PRECISION SCIENTIFIC and a complete understanding of this unit is necessary to attain these objectives.

As the ultimate user of this apparatus, it is your responsibility to-understand its proper function and operational This instruction characteristics. manual should be thoroughly read and all operators given adequate training before attempting to place this unit in Awareness of the stated service. cautions and warnings, and compliance with recommended operating parameters-together with maintenance requirementsfor safe important satisfactory operation. The unit used for should be its intended application; alterations or modifications will void the Warranty.

<u>WARNING</u>: As a routine laboratory precaution, always wear safety glasses when working with this apparatus.

This product is not intended, nor can it be used, as a sterile or patient connected device. In addition, this apparatus is not designed for use in Class I, II, or 111 locations as defined by the National Electrical Code.

Unpacking and damage

Save all packing material if apparatus is received damaged. This merchandise was carefully packed and thoroughly inspected before leaving our factory.

Responsibility for its safe delivery was assumed by the carrier upon acceptance of the shipment; therefore, claims for loss or damage sustained in transit must be made upon the carrier by the recipient as follows:

visible Loss or Damage: Note any external evidence of loss or damage on the freight bill, or express receipt, and have it signed by the carrier's agent. Failure to adequately describe such external evidence of loss or damage may result in the carrier's refusing to honor your damage claim. The form required to file such a claim will be supplied by the carrier.

Concealed Loss or Damage: Concealed loss or damage means loss or damage which does not become apparent until the merchandise has been unpacked and inspected. Should either occur, make a written request for inspection by the carrier's agent within 15 days of the delivery date; then file a claim with the carrier since the damage is the carrier's responsibility.

By following these instructions carefully, we guarantee our full support of your claim to be compensated for loss from concealed damage.

DO NOT -- FOR ANY REASON -- RETURN THIS UNIT WITHOUT FIRST OBTAINING AUTHORIZATION. In any correspondence to PRECISION SCIENTIFIC please supply the nameplate data, including catalog number and serial number.

General information

These instructions encompass the models listed below with their specific electrical characteristics:

Cat No. <u>volts</u> Hertz Amps Watts 74537 120 50/60 6.3 750 74540 220 50/60 3.4 750 74545 Gas Heat, all commercial gases

The Pensky-Martens Flash Point Tester conforms to ASTM D-93 and is designed to determine the flash point of fuel oils, lube oils, suspensions of solids, liquids that tend to form a surface film under test conditions, and other liquids.

Installation and Operation

ASTM D-93 outlines the installation and operational procedures required for the determination of flash points and should be referred to in all cases.

ASTM methods can be obtained from:

American Society for Testing Materials 1916 Race Street Philadelphia, Pennsylvania 19103

SAFETY CONSIDERATIONS AND WARNINGS:

The following guidelines are presented to supplement the existing safety rules enforced by your company:

- 1) Safety glasses should be worn by the operator and by anyone in the vicinity who could be splashed by liquid samples.
- 2) It is recommended that a fire extinguisher of Halon 1211 or CO2 (at least a 5-1b. tank size) be placed conveniently in reach of the operator to protect against fires caused by the sample which might accidentally ignite during test.

3) Service or circuit testing should be attempted only by a qualified person who has been trained with regard to the potential danger of working with live electrical circuitry.

WARNING: Disconnect the unit from the power source whenever replacing electrical components.

Electrical Connections: Important (Please Read Carefully.)

The services of a qualified technician should be used to install this unit. It should be determined that the power supply receptacle is properly polarized and grounded.

As delivered, it is supplied with a standard three-wire polarized line cord and plug for operation on 120 volts, single phase, 50/60 Hertz, or 220V, 50/60 Hertz.

WARNING: For personal safety, this unit must be properly grounded.

When a two-prong wall receptacle is encountered, it is the personal responsibility and obligation of the installer to have it replaced with a properly grounded three-prong receptacle.

WARNING: DO NOT, under any circumstances, cut or remove the third (ground) prong from the power cord. DO NOT use a two-prong adapter plug.

Determine the total amount of current presently being used by other apparatus to the circuit that will be used for this appara_us.

INSTALLATION (Contd.)

Electrical Connections: (Contd.)

It is critical that the added current demand and other equipment on the circuit not exceed the rating of the fuse or circuit breaker, in use, on this circuit.

CAUTION: Be sure that the power supply is of the same voltage as specified on the nameplate.

Explanation of controls

Electric Heat:

The Pensky-Martens Flashpoint Tester stepless transformer for variable heat control

transformer for variable heat control (from 0 to 750 watts). The reference heater dial is conveniently worked from 0 to 100. The numbers are strictly for reference to the heater wattage. To

increase the heat, turn dial counterclockwise; to decrease heat, turn dial clockwise.

Gas Heat:

Adjusting the gas heat should be done by rotating the needle valve at the base of the unit. The heat source should be centered under the opening of the heating plate.

CAUTION: Under no circumstances should products of combustion or free flame be allowed to come up around the cup.

Thermometers:

Thermometers are not supplied with the Flash Tester. When working in the 50 to 230°F (10 to 110°C) range, it is recommended that an ASTM Thermometer 9F (20 to 230°F) or 9C (-7 to 110″C) be used. When working in the 200 to 700°F (90 to 370°C) range, it Is recommended that an ASTM Thermometer 16F, (200° to 700°F) or 16C (90 to 370°C) be used.


Parts List

			T				
	Description	<u>Qty</u>	74537	74540	74545		
1	Knob, Ceramic	ob, Ceramic 2					
2	Flash Cup	1		74548			
3	Burner	1	N	/A	510056		
4	Valve, Needle	1	N.	/A	509744		
5	Orifice	1	N,	/A	509784		
6	Cover with Operating Mechanism	1		74549			
7	Tubing, Pure Gum Rubber	2 ft.		166340			
8	Screw, Orifice Pivot	1		510037			
9	Propeller, Upper (small)*	1		510032			
10	Propeller, Lower (large)*	1	ĺ	510033			
11	Shaft, Stirrer*	1		510030			
12	Orifice	1		510025			
13	Valve, Orifice	1		509738			
14	Ferrule, Thermometer	1		517407			
15	Spring*	1		510017			
16	Pulley	1		524923			
17	Flexible Shaft, Hand Stirrer	1		74553			
18	Spacer	2	1	520699			
19	Dome	1		515975			
20	Support Rod, Motor	1		507083			
21	Tube, Gas Line	1	51.	5969	520548		
22	Valve, Needle (gas line tube)	1		515966			
23	Flash Cup Holder	1		515965			
24	Thumbscrew, Flash Cup Holder	1		428915			
25	Heater, Ful Kontrol	1	536942	540094	N/A		
26	Lower Refractory with Heater*	1	61856	540095	N/A		
27	Heater, element only"	1	61876	61877	N/A		
28	Transformer (Ohmite)*	1	225086	N/A	N/A		
	Brush*	1	225087	N/A	N/A		
29	Transformer (Staco)*	1	225086	225239	N/A		
	Brush*	1	225234	225261	N/A		

^{*} Not shown.

Accessories

Slow Speed Stirrer,	115V,	60 Hz.	75765
Slow Speed Stirrer,	220V,	50 Hz.	75766



74537 & 74540

74545

Exclusive **Precision**® Warranty

PRECISION SCIENTIFIC warrants its products against defects in material or in workmanship, when used under appropriate conditions and in accordance with appropriate operating instructions for a period of no less than one (1) year from the date of delivery of the products.

Sole obligation of PRECISION SCIENTIFIC shall be to repair or replace at our option, FOB factory or locally, without charge, any part(s) that prove defective within the warranty period, provided the customer notifies PRECISION SCIENTIFIC promptly and in writing of Compensation for any such defect. labor by other than PRECISION SCIENTIFIC employees will not be our obligation. Part(s) replacement does not constitute an extension of the original warranty period.

PRECISION SCIENTIFIC MAKES NO WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR ANY OTHER WARRANTY, EXPRESSED OR IMPLIED, AS TO THE DESIGN, SALE, INSTALLATION, OR USE OF ITS PRODUCTS, AND SHALL NOT BE LIABLE FOR CONSEQUENTIAL DAMAGES RESULTING FROM THE USE OF ITS PRODUCTS.

PRECISION SCIENTIFIC will not assume responsibility for unauthorized repairs or failure as a result of unauthorized product modifications, or for repairs, replacements, or modifications negligently or otherwise improperly made or performed by persons other than PRECISION SCIENTIFIC employees or authorized representatives.

While our personnel are available to advise customers concerning general applications of all manufactured products, oral representations are not warranties with respect to particular applications and should not be relied upon if inconsistent with product specifications or the terms stated herein.

In any event, the terms and conditions contained in PRECISION SCIENTIFIC formal sales contracts shall be controlling; and any changes must be in writing and signed by an authorized executive of PRECISION SCIENTIFIC.

All defective components will be replaced without charge one (1) year from the date of delivery. There will be no charge for labor if the apparatus is returned to the factory prepaid.

Conditions and qualifications of the warranty statement shall prevail at all times.

Precision® is a registered trademark of Precision Scientific Inc.

APPENDIX A

REFERENCES

A-1 . **Scope.** This appendix contains all forms, pamphlets and technical manuals referenced in both the Air mobile and Semitrailer mounted Laboratories.

A-2. Forms.

Recommended Changes to Publications	DA Form 2028
· ·	DA Form 2028-2
Quality Deficiency Report	SF 368
Equipment Inspection and Maintenance Work Sheet	DA Form 2404
Hand Receipts	DA Form 2062
A-3. Field Manuals.	
Petroleum Testing Facilities:	
Laboratories and Kits	
inspecting and Testing Petroleum Products	
ASTM Test Method Supplement to	FM10-92C1/C2
A-4. Technical Manuals.	
Atlas-Copco Compressor	TM 10-4310-392-13&P
Alcor Jet Fuel Thermal Oxidation Tester Operating	
and Maintenance Manual	TM 106635-210-13&P
Bacharach Gas Alarm and Calibration Data	TM 106665-297-13&P
Brother Portable Typewriter	
Chemtrix Field Ph Meter	TM 10-6630-237-13&P
Elkay Manufacturing 30 GPH Cooler	TM 10-4130-240-13&P
Emcee Micro-Separometer	
Foxboro Pressure Recording Gauge	
Gammon Aqua Glo Water Detector	
Gammon Mini Monitor Fuel Sampling Kit	
Jelrus Burn-Out Furnace	
Koehler Cleveland Open Tester	
Koehler Cloud and Pour Point Chamber	
Koehler Copper Strip Corrosion Bomb Bath	
Koehler Distillation Apparatus	
Koehler Dropping Point Apparatus	
Koehler Electric Pensky-Martins Tester	
Koehler Foaming Characteristics Determination Apparatus	
Koehler Kinematic Viscosity Bath	
Koehler Tag Closed Cup Flash Tester	
Lab-Line Explosion Proof Refrigerator	
Lily Freezer	
Millipore 0M 39 Filter Holder	
Millipore Vacuum Pump	
Ohaus Harvard Trip Balance	FM 10-6670-278-13&P
Precision Gas-Oil Distillation Test Equipment	
Precision General Purpose Water Bath	IM 10-640-229-13&P

TM 10-6630-231-13&P

Precision High Temperature Bronze Block Gum Bath
Precision General Purpose Ovens
Precision Oxidation Stability Bath
Precision Pensky-Martens Flash Testers TM 10-6630-231-13&P
Precision Reid Vapor Pressure Bath
Precision W-Speed Stirrer TM 10-6640-224-13&P
Precision Universal Centrifuge TM 10-6640-230-13&P
Precision Universal Penetrometer
Sargent-Welch Vacuum Pump TM 10-4310-391-13&P
Sartorious Analytical Balance TM 10-6670-277-13&P
Scotsman Cuber
Soltec VOM-Multimeter
Teel Self-Priming Centrifugal Pump TM 10-6640-217-13&P
Teel Submersible Pump TM 10-4320-320-13&P
Texas Instrument TI-5030II Calculator TM 10-7420-210-13&P
A-5. Pamphlets.
The Army Maintenance Management System (TAMMS) DA Pam 738-750
A-6. Miscellaneous Publications.
The Army integrated Publishing and Printing Program
Clinical, College and Government Laboratories Fisher Scientific Laboratories Catalog Petroieum-Petrochemical Testing Equipment Precision Scientific Catalog

APPENDIX B

MAINTENANCE ALLOCATION CHART

Section I. INTRODUCTION

B-1. General.

- **a.** This section provides a general explanation of all maintenance and repair functions authorized at various maintenance categories.
- **b.** The Maintenance Allocation Chart (MAC) in Section II designates overall authority and responsibility for the performance of maintenance functions on the identified end item or component. The application of the maintenance functions to the end item or component will be consistent with the capacities and capabilities of the designated maintenance categories.
- **c.** Section III lists the tools and test equipment (both special tools and common tool sets) required for each maintenance function as referenced from Section II.
- **d.** Section IV contains supplemental instructions and explanatory notes for a particular maintenance function.
- B-2. Maintenance Functions. Maintenance functions will be limited to and defined as follows:
- **a.** <u>Inspect.</u> To determine the serviceability of an item by comparing its physical, mechanical, and/or electrical characteristics with established standards through examination (e.g., by sight, sound, or feel).
- **b.** <u>Test.</u> To verify serviceability by measuring the mechanical, pneumatic, hydraulic, or electrical characteristics of an item and comparing those characteristics with prescribed standards.
- **c.** <u>Service.</u> Operations required periodically to keep an item in proper operating condition, i.e., to clean (includes decontaminate, when required), to preserve, to drain, to paint, or to replenish fuel, lubricants, chemical fluids, or gases.
- **d.** <u>Adjust.</u> To maintain or regulate, within prescribed limits, by bringing into proper or exact position, or by setting the operating characteristics to specified parameters.
- **e.** <u>Align.</u> To adjust specified variable elements of an item to bring about optimum or desired performance.
- **f.** <u>Calibrate.</u> To determine and cause corrections to be made or to be adjusted on instruments or test, measuring, and diagnostic equipments used in precision measurement. Consists of comparisons of two instruments, one of which is a certified standard of knob accuracy, to detect and adjust any discrepancy in the accuracy of the instrument being compared.
- **g.** Remove/Install. To remove and install the same item when required to perform service or other maintenance functions. Install may be the act of emplacing, seating, or fixing into position a spare, repair part, or module (component or assembly) in a manner to allow the proper functioning of an equipment or system.
- **h.** <u>Replace.</u> To remove an unserviceable item and install a serviceable counterpart in its place. "Replace" is authorized by the MAC and is shown as the third position code of the SMR code.

- *i.* <u>Repair.</u> The application of maintenance services, ¹including fault location/troubleshooting, ² removal/installation, and disassembly/assembly procedures, ³ and maintenance actions ⁴ to identify troubles and restore serviceability to an item by correcting specific damage, fault, malfunction, or failure in a part, subassembly, module (component or assembly), end item, or system.
- **j.** <u>Overhaul.</u> That maintenance effort (service/action) prescribed to restore an item to a completely serviceable/operational condition as required by maintenance standards in appropriate technical publications (i.e, DMWR). Overhaul is normally the highest degree of maintenance performed by the Army. Overhaul does not normally return an item to like-new condition.
- **k.** <u>Rebuild.</u> Consists of those services/actions necessary for the restoration of unserviceable equipment to a like-new condition in accordance with original manufacturing standards. Rebuild is the highest degree of materiel maintenance applied to Army equipment. The rebuild operation includes the act of returning to zero those age measurements (hours/miles, etc.) considered in classifying Army equipment/components.

B-3. Explanation Of Columns In The MAC, Section II.

- a. <u>Column I .Group Number</u>. Column 1 lists functional group code numbers, the purpose of which is to identify maintenance significant components, assemblies, subassemblies, and modules with the next higher assembly. End item group number shall be "00."
- **b.** <u>Column 2. Component/Asse mblv</u>. Column 2 contains the names of components, assemblies, subassemblies, and modules for which maintenance is authorized.
- c. <u>Column 3. Maintenance Function.</u> Column 3 lists the functions to be performed on the item listed in column 2. (For a detailed explanation of these functions, see paragraph B-2.)
- d. <u>Column 4. Maintenance Category</u>. Column 4 specifies, by the listing of a work time figure in the appropriate subcolumn(s), the category of maintenance authorized to perform the function listed in column 3. This figure represents the active time required to perform that maintenance function at the indicated category of maintenance. If the number or complexity of the tasks within the listed maintenance function vary at different maintenance categories, appropriate work time figures will be shown for each category. The work time figure represents the average time required to restore an item (assembly, subassembly, component, module, end item, or system) to a serviceable condition under typical field operating conditions. This time includes preparation time (including any necessary disassembly/ assembly time), troubleshooting/fault location time, and quality assurance/quality control time in addition to the time required to perform the specific tasks identified for the maintenance functions authorized in the maintenance allocation chart. The symbol designations for the various maintenance categories are as follows:

¹ Services - inspect, test, service, adjust, align, calibrate, and/or replace.

² Fault locate/troubleshoot -the process of investigating and detecting the cause of equipment malfunctioning; the act of isolating a fault within a system or unit under test (UUT).

³ Disassemble/assemble- encompasses the step-by-step taking apart (or breakdown) of a spare/functional group coded item to the level of its least componency identified as maintenance significant (i. e., assigned an SMR code) for the category of maintenance under consideration.

⁴ Actions - welding, grinding, riveting, straightening, facing, remachining, and/or resurfacing.

C	Operator/Crew
0	Unit Maintenance
F	Direct Support Maintenance
H	General Support Maintenance
D	Depot Maintenance

- e. <u>Column 5. Tools and Equipment.</u> Column 5 specifies, by code, those common tool sets (not individual tools) and special tools," TM DE, and support equipment required to perform the designated function.
- f. <u>Column 6. Remarks.</u> This column shall, when applicable, contain a letter code, in alphabetic order, which shall be keyed to the remarks contained in section IV.
- B-4. Explanation Of Columns In Tool And Test Equipment Requirements, Section III.
- a. <u>Column I. Reference Code.</u> The tool and test equipment reference code correlates with a code used in the MAC, section II, column 5.
- b. <u>Column 2. Maintenance Category.</u> The lowest category of maintenance authorized to use the tool or test equipment.
 - c. <u>Column 3. Nomenclature.</u> Name or identification of the tool or test equipment.
 - d. Column 4. National Stock Number. The National stock number of the tool or test equipment.
 - e. <u>Column 5. Tool Number</u>. The manufacturer's part number.
- B-5. Explanation Of Columns In Remarks, Section IV.
 - a. <u>Column I. Reference Code</u>. The code recorded in column 6, Section II.
- b. <u>Column 2. Remarks.</u> This column lists information pertinent to the maintenance function being performed as indicated in the MAC, section II.

Section II. MAINTENANCE ALLOCATION CHART

(1) GROUP NUMBER	(2) COMPONENT/ ASSEMBLY	(3) MAINTENANCE FUNCTION		MAINT NIT	(4 ENAN DS F	ICE L	EVEL DEPOT D	(5) TOOLS AND EQUIPMENT	(6) REMARKS
	TESTER, FLASH POINT (CLOSED)	INSPECT REPLACE REPAIR	0.3	0.3 1.5	_			1, 2	A

TM10-6630-231-13&P

SECTION III. TOOL AND TEST EQUIPMENT REQUIREMENTS

FOR

MAINTENANCE ALLOCATION CHART

(1) (2) (3) (4) (5) TOOL/TEST

EQUIP. MAINTENANCE NOMENCLATURE TOOL NSN REF CODE CATEGORY NUMBER

TOOL KIT, GENERAL AUTOMOTIVE 5180-00-177-7033 (50980)

SC 5180-90-CL-N26

MULTIMETER, 0-500V 6625-00-691-2453 2

SECTION IV. REMARKS

REFERENCE

CODE REMARKS

REPLACE DEFECTIVE PARTS

APPENDIX C

COMPONENTS OF END ITEM AND BASIC ISSUE ITEMS LISTS Section I. INTRODUCTION

C-1. Scope.

This appendix lists components of end item and basic issue items for the Electric Pensky-Martens Tester to help you inventory items required for safe and efficient operation.

C-2. General.

The Components of End Item and Basic Issue Items Lists are divided into the following sections:

- a. <u>Section II. Components of End Item</u>. This listing is for informational purposes only, and is not authority to requisition replacements. These items are part of the end item, but are removed and separately packaged for transportation or shipment. As part of the end item, these items must be with the end item whenever it is issued or transferred between property accounts. Illustrations are furnished to assist you in identifying the items.
- b. <u>Section III. Basic Issue Items</u>. These are the minimum essential items required to place the Electric Pensky-Martens Tester in operation, to operate it, and to perform emergency repairs. Although shipped separately packaged, BII must be with the shelter during operation and whenever it is transferred between property accounts. The illustrations will assist you with hard-to-identify items. This manual is your authority to reques/requisition replacement BII, based on TOE/MTOE authorization of the end item.

C-3. Explanation of Columns.

The following provides an explanation of columns found in the tabular listings:

- **a.** <u>Column (1) Illustration Number (Illus Number)</u>. This column indicates the number of the illustration in which the item is shown.
- **b.** <u>Column (2) National Stock Number.</u> Indicates the National stock number assigned to the item and will be used for requisitioning purposes.
- c. <u>Column (3) Description</u>. Indicates the Federal item name and, if required, a minimum description to identify and locate the item. The last line for each item indicates the CAGEC (in parentheses) followed by the part number.
- d. <u>Column (4) Unit of Measure (U/M)</u>. Indicates the measure used in performing the actual operational/maintenance function. This measure is expressed by a two-character alphabetical abbreviation (e.g., ea, in, pr).
- **e.** <u>Column (5) Quantity required (QTY RQR)</u>. Indicates the quantity of the item authorized to be used with/on the equipment.

TM10-6630-321-13&P

SECTION II. COMPONENTS OF END ITEM

(1) (2) (3) (4) (5)

NATIONAL STOCK DESCRIPTION USABLE

ILLUS NUMBER CAGEC AND PART NUMBER ON CODE U/M QTY

REGULATOR EA 2

(05083) 0023-1492

SECTION III. BASIC ISSUE ITEMS NOT APPLICABLE

APPENDIX D ADDITIONAL AUTHORIZATION LIST

NOT APPLICABLE

APPENDIX E

EXPENDABLE/DURABLE SUPPLIES AND MATERIALS LIST

Section I. INTRODUCTION

E-1. **Scope.** This listing is for informational purposes only and is not authority to requisition the listed items. These items are authorized to you by CTA 50-970, Expendable/Durable Items (except medical, class V, repair parts, and heraldic items).

E-2. Explanation of Columns.

- a. <u>Column (1) Item Number</u>. This number is assigned to the entry in the listing and is referenced in the narrative instructions to identify the material (e.g., Use cleaning compound, item 5, appendix C).
- b. <u>Column (2) Level.</u> This column identifies the lowest level of maintenance that requires the listed item.
 - C Operator/Crew
 - O Unit Maintenance
 - F Direct Support Maintenance
 - H General Support Maintenance
- c. <u>Column [3] National Stock Number.</u> This is the National stock number assigned to the item; use it to request or requisition the item.
- d. <u>Column (4) Description</u>. Indicates the Federal item name, and, if required, a description to identify the item. The last line for each item indicates the Commercial and Government Entity Code (CAGEC) in parentheses followed by the part number.
- e. <u>Column (5) Unit of Measure (U/M)</u>. Indicates the measure used in performing the actual maintenance function. This measure is expressed by a two-character alphabetical abbreviation (e.g., EA,IN, PR). if the unit of measure differs from the unit of issue, requisition the lowest unit of issue that will satisfy your requirements.

Section II. EXPENDABLE/DURABLE SUPPLIES AND MATERIALS LIST

(1) Item Number	(2) Level	(3) National Stock Number	(4) Description	(5) U/M
	С	6830-00-269-4299	PROPANE: CYLINDER DISPOSABLE, ODORIZED (SMALL) 1 LB. CYL (80244) BB-6-40, TYPE III	EA

By Order of the Secretary of the Army:

CARL E. VUONO

General, United States Army Chief of Staff

Official:

THOMAS F. SIKORA

Brigadier General, United States Army The Adjutant General

DISTRIBUTION:

To be distributed in accordance with DA Form 12-21A, Operator, Unit and Direct Support Maintenance requirements for Laboratory, Air Mobile, Aviation Fuel and Laboratory, Petroleum, MTD

☆ U.S. GOVERNMENT PRINTING OFFICE: 1991 554-123/20088



SOMETHING WRONG WITH THIS PUBLICATION?

THEN.. JOT DOWN THE DOPE ABOUT IT ON THIS FORM, CAREFULLY TEAR IT OUT, FOLD IT AND DROP IT IN THE MAIL!

FROM: (PRINT YOUR UNIT'S COMPLETE ADDRESS)

PFC JOHN DOE

COA, 34 ENGINEER BN

FT. LEONARDWOOD, MO 63108

DATE SENT

PUBLICATION NUMBER

TEAR ALONG PERFORATED LINE

TM 10-6630-231-13&P

PUBLICATION DATE 28 Sep 1990

PUBLICATION TITLE Koehler Electric Pensky-Martens Tester

	_				
BE EXACTPIN-POINT WHERE IT IS		IN THIS	S SPACE TELL WHAT I	s wrong	
PAGE PAR NO. GRA		TABLE NO.		HAT SHOULD BE DON	
6 2			Ju	line 6g	paragraph 2-10 the
	٦		m	ranual s	totes the engine has
			6	Cylendus.	totes the engine has
			se	tonly ho	s 4 Cylenders.
			ck	enge the .	manual to show L
	j			fenders.	
			0		
81	4-3		Ca	llaut 16 o	n figure 4-3 is
			0	interior of	a pour, money
		i I	15	lieute 4	- 3, item 16 is celled
			a	shim -	Please Correct
			01	re to the	Other.
125		20	9	redered o	gasket, item
'45 ^			19	on Sign	ne B-16 leg NSN
				910-00-76	2-3001. I get a
				abot lo	t it dream it fit.
			g	S A	and what
			X	upply X	eyo I got what
			A	ordered,	so the NSN is
			W	rong. The	ease give me a
PRINTED NAME. G	RADE OR TITLE	AND TELEPH	IONE NUME	SIGN HE	RE ONE DA

PRINTED NAME, GRADE OR TITLE, AND TELEPHONE NUMBER

JOHN DOE, PFC (268) 317.7111

JOHN DOE

DA 1 JUL 79 2028-2

PREVIOUS EDITIONS ARE OBSOLETE.

DRSTS-M Overprint 1, 1 Nov 80

P.S.--IF YOUR OUTFIT WANTS TO KNOW ABOUT YOUR RECOMMENDATION MAKE A CARBON COPY OF THIS AND GIVE IT TO YOUR HEADQUARTERS.

TEAR ALONG PERFORATED LINE



FOLD BACK

DEPARTMENT OF THE ARMY

OFFICIAL BUSINESS

COMMANDER
U.S. ARMY TROOP SUPPORT COMMAND
ATTN: AMSTR-MCTS
4300 GOODFELLOW BOULEVARD
ST. LOUIS, MO 63120-1798

RECOMMENDED CHANGES TO EQUIPMENT TECHNICAL PUBLICATIONS



SOMETHING WRONG WITH THIS PUBLICATION?

THEN. JOT DOWN THE DOPE ABOUT IT ON THIS FORM, CAREFULLY TEAR IT OUT. FOLD IT AND DROP IT IN THE MAIL! FROM: (PRINT YOUR UNIT'S COMPLETE ADDRESS)

DATE SENT

PUBLICATION NUMBER

ALONG PERFORATED UNF

PUBLICATION DATE

PUBLICATION TITLE Koehler Electric Pensky-Martens Tester

TH 10-6630-231-13&P

28 Sep 1990

PIN-POINT WHERE IT IS BE EXACT. IN THIS SPACE TELL WHAT IS WRONG PARA-GRAPH FIGURE TABLE AND WHAT SHOULD BE DONE ABOUT IT:

SIGN HERE.

DA 1 JUL 79 2028-2

PRINTED NAME, GRADE OR TITLE, AND TELEPHONE NUMBER

PREVIOUS EDITIONS ARE OBSOLETE. DRSTS-M Overprint 2, 1 Nev 80

P.S.--IF YOUR OUTFIT WANTS TO KNOW ABOUT YOUR RECOMMENDATION MAKE A CARBON COPY OF THIS AND GIVE IT TO YOUR HEADQUARTERS.

REVERSE OF DA FORM 2028-2 REVERSE OF DRSTS-M OVERPRINT 2, 1 NOV 80

TEAR ALONG PERFORATED LINE



FOLD BACK

DEPARTMENT OF THE ARMY

OFFICIAL BUSINESS

COMMANDER
U.S. ARMY TROOP SUPPORT COMMAND
ATTN: AMSTR-MCTS
4300 GOODFELLOW BOULEVARD
ST. LOUIS, MO 63120-1798

RECOMMENDED CHANGES TO EQUIPMENT TECHNICAL PUBLICATIONS



SOMETHING WRONG WITH THIS PUBLICATION?

THEN. JOT DOWN THE DOPE ABOUT IT ON THIS FORM, CAREFULLY TEAR IT OUT, FOLD IT AND DROP IT IN THE MAIL'

FROM: (PRINT YOUR UNIT'S COMPLETE ADDRESS)

DATE SENT

PUBLICATION NUMBER

ALONG PERFORATED LINE

PUBLICATION DATE

PUBLICATION TITLE Koehler Electric Pensky-Martens Tester

Tri 10-6630-231-13&P

BE EXACT. PIN-POINT WHERE IT IS

28 Sep 1990

IN THIS SPACE TELL WHAT IS WRONG AND WHAT SHOULD BE DONE ABOUT IT:

FIGURE NO TABLE PAGE NO PARA-GRAPH PRINTED NAME, GRADE OR TITLE, AND TELEPHONE NUMBER SIGN HERE

DA 1 JUL 79 2028-2

PREVIOUS EDITIONS
ARE OBSOLETE.

DRSTS-M Overprint 2, 1 Nev 80.

P.S.--IF YOUR OUTFIT WANTS TO KNOW ABOUT YOUR RECOMMENDATION MAKE A CARBON COPY OF THIS AND GIVE IT TO YOUR HEADQUARTERS.

REVERSE OF DA FORM 2028-2 REVERSE OF DRSTS-M OVERPRINT 2, 1 NOV 80

TEAR ALONG PERFORATED LINE



FOLD BACK

DEPARTMENT OF THE ARMY

OFFICIAL BUSINESS

COMMANDER U.S. ARMY TROOP SUPPORT COMMAND ATTN: AMSTR-MCTS 4300 GOODFELLOW BOULEVARD ST. LOUIS, MO 63120-1798

The Metric System and Equivalents

Linear Measure

1 centimeter = 10 millimeters = .39 inch 1 decimeter = 10 centimeters = 3.94 inches 1 meter = 10 decimeters = 39.37 inches 1 dekameter = 10 meters = 32.8 feet 1 hectometer = 10 dekameters = 328.08 feet 1 kilometer = 10 hectometers = 3,280.8 feet

Weights

1 centigram = 10 milligrams = .15 grain 1 decigram = 10 centigrams = 1.54 grains 1 gram = 10 decigram = .035 ounce 1 dekagram = 10 grams = .35 ounce 1 hectogram = 10 dekagrams = 3.52 ounces 1 kilogram = 10 hectograms = 2.2 pounds 1 quintal = 100 kilograms = 220.46 pounds 1 metric ton = 10 quintals = 1.1 short tons

Liquid Measure

1 centiliter = 10 milliters = .34 fl. ounce 1 deciliter = 10 centiliters = 3.38 fl. ounces 1 liter = 10 deciliters = 33.81 fl. ounces 1 dekaliter = 10 liters = 2.64 gallons 1 hectoliter = 10 dekaliters = 26.42 gallons 1 kiloliter = 10 hectoliters = 264.18 gallons

Square Measure

1 sq. centimeter = 100 sq. millimeters = .155 sq. inch 1 sq. decimeter = 100 sq. centimeters = 15.5 sq. inches 1 sq. meter (centare) = 100 sq. decimeters = 10.76 sq. feet 1 sq. dekameter (are) = 100 sq. meters = 1,076.4 sq. feet 1 sq. hectometer (hectare) = 100 sq. dekameters = 2.47 acres 1 sq. kilometer = 100 sq. hectometers = .386 sq. mile

Cubic Measure

1 cu. centimeter = 1000 cu. millimeters = .06 cu. inch 1 cu. decimeter = 1000 cu. centimeters = 61.02 cu. inches 1 cu. meter = 1000 cu. decimeters = 35.31 cu. feet

Approximate Conversion Factors

To change	To	Multiply by	To change	To	Multiply by
inches	centimeters	2.540	ounce-inches	newton-meters	.007062
feet	meters	.305	centimeters	inches	.394
yards	meters	.914	meters	feet	3.280
miles	kilometers	1.609	meters	yards	1.094
square inches	square centimeters	6.451	kilometers	miles	.621
square feet	square meters	.093	square centimeters	square inches	.155
square yards	square meters	.836	square meters	square feet	10.764
square miles	square kilometers	2.590	square meters	square yards	1.496
acres	square hectometers	.405	square kilometers	square miles	.386
cubic feet	cubic meters	.028	square hectometers	acres	2.471
cubic yards	cubic meters	.765	cubic meters	cubic feet	35.315
fluid ounces	milliliters	29,573	cubic meters	cubic yards	1.308
pints	liters	.473	milliliters	fluid ounces	.034
quarts	liters	.946	liters	pints	2.113
gallons	liters	3.785	liters	quarts	1.057
ounces	grams	28.349	liters	gallons	.264
pounds	kilograms	.454	grams	ounces	.035
short tons	metric tons	.907	kilograms	pounds	2.205
pound-feet	néwton-meters	1.356	metric tons	short tons	1.102
pound-inches	newton-meters	.11296			•

Temperature (Exact)

٥F	Fahrenheit	5/9 (after	Celsius	$^{\circ}\mathbf{C}$
	temperature	subtracting 32)	temperature	

PIN: 046473-000

This fine document...

Was brought to you by me:



<u>Liberated Manuals -- free army and government manuals</u>

Why do I do it? I am tired of sleazy CD-ROM sellers, who take publicly available information, slap "watermarks" and other junk on it, and sell it. Those masters of search engine manipulation make sure that their sites that sell free information, come up first in search engines. They did not create it... They did not even scan it... Why should they get your money? Why are not letting you give those free manuals to your friends?

I am setting this document FREE. This document was made by the US Government and is NOT protected by Copyright. Feel free to share, republish, sell and so on.

I am not asking you for donations, fees or handouts. If you can, please provide a link to liberatedmanuals.com, so that free manuals come up first in search engines:

Free Military and Government Manuals

- SincerelyIgor Chudovhttp://igor.chudov.com/
- Chicago Machinery Movers